


## TYPE EVD-VE110



**Pumping speed curve**



Technical Data		TYPE EVD-VE110	
		50 Hz	60 Hz
Nominal displacement	m <sup>3</sup> /h	1,7	2,0
Ultimate pressure	mbar	0,05	0,05
Motor power	W	185	185
Oil capacity	ltr	0,3	0,3
Weight	kg	7	7
Admissible ambient temp.	°C	10 to 40	10 to 40
Connections	Flare	¼"	¼"
Dimensions LxWxH	mm	285x100x230	285x100x230

Ordering Information	EVD-VE110
Single phase motor 220-240V, 50/60 Hz (IEC)	447000
Dual Voltage motor 220V, 50 Hz & 110V, 60 Hz	447100
Exhaust filter	410900
Accessories	See page 13