


TYPE EVM-D1.2



Pumpdown time for 20 ltr vessel



Technical Data		EVM-D1.2	EVM-D1.2C
		Standard version	Chemical version
Nominal displacement	m ³ /h	1,2	
Ultimate vacuum	mbar	50	
Motor power	W	160	
Material Diaphragm & Valves		NBR	PTFE
Material Pump Head coating		none	PTFE
Weight	kg	10	
Noise level	dB	< 60	
Admissible ambient temp.	°C	7 to 40	
Connections	mm	6	
Dimensions LxWxH	mm	300x120x235	

Ordering Information	EVM-D1.2	EVM-D1.2C
Single phase motor 220-240V, 50 Hz (IEC)	623000	623100
Accessories	See page 9	