


TYPE EVM-S1.2



Pumpdown time for 10 ltr vessel



Technical Data		EVM-S1.2	EVM-S1.2C
		Standard version	Chemical version
Nominal displacement	m ³ /h	1,2	
Ultimate vacuum	mbar	200	
Motor power	W	160	
Material Diaphragm & Valves		NBR	PTFE
Material Pump Head coating		none	PTFE
Weight	kg	7,5	
Noise level	dB	< 60	
Admissible ambient temp.	°C	7 to 40	
Connections	mm	6	
Dimensions LxWxH	mm	215x120x235	

Ordering Information	EVM-S1.2	EVM-S1.2C
Single phase motor 220-240V, 50 Hz (IEC)	622000	622100
Accessories	See page 9	